

# HydraMed Forte®

**Moisturising, lubricant sterile ophthalmic solution**  
**Sodium Hyaluronate 0.4% and TS-Polysaccharide 0.2%**  
 30 x 0.5 ml daily use preservative free reclosable vials

**Composition:** Sodium Hyaluronate (0.4% w/v), TS-Polysaccharide (0.2% w/v), Mannitol, Sodium citrate, Citric acid monohydrate, Water for injections.

#### What is it?

**HydraMed Forte** is a sterile ophthalmic solution containing Sodium Hyaluronate and TS-Polysaccharide.

**Phosphate and preservative free.**

#### How is it presented?

**HydraMed Forte** is presented as an eye drops solution in 30 daily use re-closable vials of 0.5 ml.

#### When should it be used?

**HydraMed Forte** is an ophthalmic solution of TS-Polysaccharide and high-molecular weight Sodium Hyaluronate, produced through biotechnological synthesis, used to obtain long lasting relief from dry eyes.

Sodium Hyaluronate is a substance that originates mucomimetic solutions, with optimal lubricant properties due to the viscosity of the solutions that result in a moisturising effect and stabilization of the tear film.

TS-Polysaccharide is a natural substance extracted from the seeds of the Tamarindus indica plant, useful in the protection of the corneo-conjunctival surface. TS-Polysaccharide forms a protective film on the ocular surface that helps to recover the physiologic balance of the tear film and the ocular surface. The combination of these two substances generates a synergic action on the corneo-conjunctival surface that is protected from irritation due to environmental elements (wind, sun, smoke, dry air), visual stress (excessive light and prolonged use of computers) or mechanical factors (contact lenses use, after ocular surgery), providing a long lasting relief.

The concentration of the mixture is calibrated in such a way as to ensure the lubrication even in cases of more intense discomfort.

#### How and how much should it be used?

- Open the aluminum bag and separate one vial from the strip of vials;
- Open it by turning and pulling the cap;
- Put one drop in the eye by squeezing gently the body of the vial, one or more times daily, according to necessity.

HydraMed Forte is filled in re-closable vials and the remaining product can be used up to 12 hours.


HydraMed Forte is compatible with contact lenses and it is not necessary to remove them before the ophthalmic solution is applied.

#### Precautions and warnings

- The product is for external ophthalmic use only.
- Do not use if the container is damaged.
- Do not use in case of known incompatibility to any of the ingredients of the preparation.
- Do not use after the expiry date shown on the pack.
- Do not touch the eye with the tip of the mini-container.
- After administration, in rare cases, a slight temporary blurring of vision may be observed due to the viscosity of the solution; it is therefore advisable to wait until the blurring has disappeared before undertaking any activity that requires clear vision.
- Leave at least 10-15 minutes before administration of any other ophthalmic drugs.
- Once opened, each vial can be used for up to 12 hours.
- Keep out of the sight and reach of children.
- Store below 25°C.
- Do not swallow.

**Manufacturer:** Farmigea S.p.A. - Via G.B. Oliva, 6/8 - 56121 Pisa, Italy

**Date of last revision of the text:** March 2019



CE 0546

STERILE A

Dispose of properly after use.

